

Common Name	Internal Diameter
Capillary LC	300µm - 1.0mm
Micro-bore	1.0mm
Narrow-bore	2.1mm
Mid-bore/Solvent Saver	3.0mm
Analytical Bore	4.6mm
Semi-Prep	7.8-21.2mm
Preparative	up to 3"
Process Scale	over 3" diameters



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com