

## The common names of HPLC column diameters - Primer

<u>Common Name</u>	<u>Internal Diameter</u>
Capillary LC	300µm - 1.0mm
Micro-bore	1.0mm
Narrow-bore	2.1mm
Mid-bore/Solvent Saver	3.0mm
Analytical Bore	4.6mm
Semi-Prep	7.8-21.2mm
Preparative	up to 3"
Process Scale	over 3" diameters



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 09-16-2024