



The common names of HPLC column diameters – Primer

Common Name	Internal Diameter
Capillary LC	300µm – 1.0mm
Micro-bore	1.0mm
Narrow-bore	2.1mm
Mid-bore/Solvent Saver	3.0mm
Analytical Bore	4.6mm
Semi-Prep	7.8-21.2mm
Preparative	up to 3"
Process Scale	over 3" diameters



MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com