



Can the all silicone cap mats be sterilized – FAQ

Yes, the silicone cap mats can be sterilized using heat at 121°C for 15-30 minutes at 16-17 psi.

Click [HERE](#) for Silicone rubber sealing mat ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 03-05-2024