



## Which ferrules are best when using soft polymeric tubing - FAQ

Kel-F™ (*PCTFE which is similar to PTFE*) and Tefzel™ (*ETFE*) ferrules give the best results when using PTFE, FEP, PFA or ETFE tubing.

[More information about PEEK tubing](#)

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024