



## Is PEEK tubing or PEEK fittings autoclavable – FAQ

Yes, many scientists have successfully used PEEK tubing and or fittings after autoclaving with gamma irradiation, heat or ethylene oxide sterilization.

Click [HERE](#) for PEEK Tubing Information

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-15-2024