



Is PEEK tubing or PEEK fittings autoclavable - FAQ

Yes, many scientists have successfully used PEEK tubing and or fittings after autoclaving with gamma irradiation, heat or ethylene oxide sterilization.

Click [HERE](#) for PEEK Tubing Information

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024