



Are stainless steel tubing or fittings autoclavable – FAQ

Yes, many scientists have successfully used 316 stainless steel after autoclaving tubing and or fittings with gamma irradiation, heat or ethylene oxide sterilization methods.

Click [HERE](#) for stainless steel tubing ordering information

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com