



## Are stainless steel tubing or fittings autoclavable – FAQ

Yes, many scientists have successfully used 316 stainless steel after autoclaving tubing and or fittings with gamma irradiation, heat or ethylene oxide sterilization methods.

Click [HERE](#) for stainless steel tubing ordering information

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)