



What is the pH range of 316 stainless steel laboratory tubing fittings and sparging stones - FAQ

The useful pH range of 316 stainless steel tubing, fittings and sparging stones is 1 - 14.

Click [HERE](#) for more information about stainless steel tubing

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024