



What is the pH range of 316 stainless steel laboratory tubing fittings and sparging stones - FAQ

The useful pH range of 316 stainless steel tubing, fittings and sparging stones is 1 - 14.

Click [HERE](#) for more information about stainless steel tubing

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com