



## What is the useful temperature range of 316 stainless steel laboratory products - FAQ

The useful temperature range of stainless steel (SS) related products such as column hardware, tubing and HPLC fittings is up to 240°C.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 03-05-2024