



## Are RSA autosampler vials pyrogen free - FAQ

YES. The RSA™ or Reduced Surface **Activity** autosampler vials are pyrogen free when manufactured, as they are produced sterile and are never touched by “human hands”. They are then packaged, including shrink wrapped in a sterile, class 10,000 clean room environment.

We do not classify these vials as sterile for sale since the packaging is not sterilized before use. However, this can give you a very good idea of the cleanliness if you wished to sterilize them.

The entire manufacturing process is ISO Certified.

[More information about RSA vials, inserts and caps](#)

