



## Are the Cogent TYPE-C Silica based stationary phases available for sale as bulk materials – FAQ

No, the Cogent 4um and 2.0 TYPE-C Silica™ based stationary phases are not available for sale as unpacked bulk material.

TYPE-C Silica™ is a very unique material and to properly pack columns or capillaries requires proprietary technology. For this reason, it is only sold in a pre-packed column format. However, 10um particles can be made available upon request.



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 04-05-2024