

## Are 10um pore size available in the Cogent Column Filters – FAQ

**No, we only offer the pre-column, Cogent Column Filters in the following pore sizes:**

- 2um stainless steel mesh: Catalog Number [49650-020](#)
- 0.5um stainless steel mesh: Catalog Number [49650-005](#)

*HPLC Column Protection*



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)