



## Are AQ brand nylon syringe filters compatible with DMSO – FAQ

**Yes, MicroSolv's AQ™ brand of syringe filters for dissolution testing and HPLC are made with a polypropylene housing which is DMSO compatible.** Our nylon membranes, both 0.45um and 0.22um are fully compatible with solutions containing DMSO. Nylon is also compatible with esters, alcohols, weak acids and bases.

Click [HERE](#) for ordering information and pictures of AQ™ brand of syringe filters.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)