

## What is the structure of the Cogent Phenyl Hydride bonded ligand – FAQ

**The Cogent Phenyl Hydride phase is a 4 carbon chain with a phenyl group at the end attached directly to the silica hydride surface with a single, direct silicon-carbon bond.**

With bonds in traditional silica B based materials, the organo-silane starting material is Si-O-Si-C4-Phenyl and this reagent is reacted with the silica to form the final stationary phase which includes a Si-O-Si-C4-Phenyl ligand.

This is in contrast with the TYPE-C phases in that the Si-C bond is formed directly as part of the reaction resulting in Si-C4-Phenyl without the oxygen in the link which is very susceptible to hydrolysis.

[More information on Phenyl Hydride HPLC Columns](#)



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 06-05-2024