

## Which vials compatible with SPME – FAQ

**The following two items are compatible with Solid Phase Micro-Extraction (SPME):**

Catalog # [95020-45S-T](#): Magnetic Screw Cap for Headspace. Silver, 8mm Center Hole, Silicone Transparent Blue/PTFE White, 45 Shore A, 1.3mm Thick. **Precision** Threads. 100/PK.

Catalog # [95010-0CV](#): Vials, head space, 10ml, 46 x 22.5mm, screw top, clear borosilicate glass, round bottom. CTC Compliant. **Precision** threads 100/PK.

However, the following caps are especially suited for SPME work: Catalog # [95020-21S-T](#): Caps, magnetic screw, for Head Space. Silver 8mm center hole with a pre slit (star shape) white silicone rubber / red PTFE septa, 55 shore A, 1.5mm thick. 18mm **precision** threads. 100/PK.

*Compatibility with specific manufacturers' autosampler instrumentation is as follows:*

*Agilent G1888A*

*Agilent/CTC Leap Combi Pal*

*CTC Combi Pal HS Mode*

*Gerstel/Varian MPS 2*

*CTC Combi Pal SPME Mode*

*Turbomax HS16*

*Turbomax HS40/40XI*

*Turbomax HS40/Trap*

*Turbomax HS110 HS110/Trap*

*Shimadzu AOC 5000 HS*

*Shimadzu AOC 5000 SPME*

[Ordering information for Headspace Vials and Caps](#)

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024