

The following two items are compatible with Solid Phase Micro-Extraction (SPME):

Catalog # [95020-45S-T](#): Magnetic Screw Cap for Headspace. Silver, 8mm Center Hole, Silicone Transparent Blue/PTFE White, 45 Shore A, 1.3mm Thick. **Precision** Threads. 100/PK.

Catalog # [95010-0CV](#): Vials, head space, 10ml, 46 x 22.5mm, screw top, clear borosilicate glass, round bottom. CTC Compliant. **Precision** threads 100/PK.

However, the following caps are especially suited for SPME work: Catalog # [95020-21S-T](#): Caps, magnetic screw, for Head Space. Silver 8mm center hole with a pre slit (star shape) white silicone rubber / red PTFE septa, 55 shore A, 1.5mm thick. 18mm **precision** threads. 100/PK.

Compatibility with specific manufacturers' autosampler instrumentation is as follows:

Agilent G1888A

Agilent/CTC Leap Combi Pal

CTC Combi Pal HS Mode

Gerstel/Varian MPS 2

CTC Combi Pal SPME Mode

Turbomax HS16

Turbomax HS40/40XI

Turbomax HS40/Trap

Turbomax HS110 HS110/Trap

Shimadzu AOC 5000 HS

Shimadzu AOC 5000 SPME

Ordering information for Headspace Vials and Caps

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com