



## For the bonded or coated CE capillaries do you have other lengths available - FAQ

**The bonded or coated capillaries are only available in 1M lengths.**

The reason is that it is very difficult to push the bonded reagent through a 50um or 75um tube and get proper coating. Click [HERE](#) for ordering information Controlled Flow and Zero EOF Capillaries

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024