

According to the WikiAnswers website, an “**assay**” method provides an experimentally determined value for the content or potency of an analyte in a sample.

In contrast, a “**purity**” methods provides an accurate quantitative statement of all impurities in a pharmaceutical sample.

This includes related substances, residual solvents, etc. Please refer to the International Conference on Harmonization (ICH) guidelines for more detail.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com