

The difference between potency and purity in drug products – Tips & Suggestions

Please refer to the International Conference on Harmonization (ICH) guidelines for more detail.

A quick definition would be:

–**Potency** is a measure of drug activity expressed in terms of the amount required to produce an effect of given intensity.

–**Purity** is a measure of the amount of API present in a sample compared to those of related substances, impurities, residual solvents, etc.



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-08-2024