



What is the inner diameter tolerance of the 50um ID for CE capillaries – FAQ

The tolerance of the inner diameter (ID) is $\pm 3 \mu\text{m}$ for all CE capillaries sold by MicroSolv.

[CE Capillary Product Page](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-26-2024