


What are the size and materials of the mobile phase delivery bottle caps and O-rings - FAQ

The mobile phase delivery bottle caps are made of polypropylene with a PTFE insert. The O-rings are made of Viton® and coated with PTFE. The caps, inserts, and O-rings are made for use with GL 45 screw thread bottles. The ring has an **outer diameter** of approximately **36.4 mm** and an **inner diameter** of approximately **28.6 mm**.



OD: 36.4mm
ID: 28.6 mm

Click [HERE](#) for mobile phase delivery bottle caps ordering information and pictures.