



## What are the size and materials of the mobile phase delivery bottle caps and O-rings – FAQ

The mobile phase delivery bottle caps are made of polypropylene with a PTFE insert. The O-rings are made of Viton® and coated with PTFE. The caps, inserts, and O-rings are made for use with GL 45 screw thread bottles. The ring has an **outer diameter** of approximately **36.4 mm** and an **inner diameter** of approximately **28.6 mm**.



Click [HERE](#) for mobile phase delivery bottle caps ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-10-2024