


What are the size and materials of the mobile phase delivery bottle caps and O-rings - FAQ

The mobile phase delivery bottle caps are made of polypropylene with a PTFE insert. The O-rings are made of Viton® and coated with PTFE. The caps, inserts, and O-rings are made for use with GL 45 screw thread bottles. The ring has an **outer diameter** of approximately **36.4 mm** and an **inner diameter** of approximately **28.6 mm**.



OD: 36.4mm
ID: 28.6 mm

Click [HERE](#) for mobile phase delivery bottle caps ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com