

Calculate the internal volume of tubing – How to

To calculate the internal volume (V), you need to know the length of tubing (L) and the inner diameter (ID). You then use the formula for volume of a cylinder.

Convert L and ID into cm first. (1 inch = 2.54 cm). This gives V in cm³. 1 cm³ = 1 mL. You can then convert into µL if desired by 1000µL = 1mL.

$$V = \pi (ID/2)^2 L$$



[Tubing Product Page](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-16-2024