

## What is the largest ID for PEEK columns that you offer - FAQ

**The largest I.D. available is 7.5mm.**

The corresponding lengths offered for a column of this ID are 30, 50, 100, 150, 250, 300 mm. All of these columns are 5/8" O.D.



*[PEEK Hardware Product Page](#)*

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-03-2024