



## What is the difference between ferrules made of stainless steel vs stainless steel and PEEK hybrids – FAQ

### **Stainless steel ferrules are used exclusively when:**

- Absolute maximum operating pressure is required, in excess of 15,000 PSI.
- The connection is 'permanent'; that is, the connection will not be made and broken often, such as when changing columns.

### **The Stainless Steel / PEEK hybrid ferrule (*called a high performance ferrule*) is:**

- Typically limited to 15,000 PSI with multi cone ferrule fittings
- Limited to 10,000 PSI with a single cone HPLC nut
- Can be re-used a limited number of times if not over torqued, or if used on SS tubing

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-05-2024