

What is included with the Whatman brand SAX HPLC column kit 63185-SCK - FAQ

The product consists of:

1. 1 each Partisil™ 10um SAX 250 x 4.6mm id standard analytical HPLC column
2. 1 each Solvecon™ pre-column kit for above: 15g Solvecon™ silica, empty 250 x 4.6mm id column (including frits), 2 each nuts & ferrules

This pre-column is designed to be a solvent conditioning column, like the old saturation columns that go in line between the pump and the Injector. It helps to saturate the mobile phase with silica and minimize analytical column dissolution problems for the Whatman SAX columns.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com