



Are the BASIK brand headspace vials free of endotoxins - FAQ

These vials are not tested for the presence of endotoxins but since they are manufactured at 600°C, bacteria or organic matter will not likely exist after manufacturing.

While the vials are packed directly after cooling, they are *not* certified to be “untouched by human hands” and may or may not contain endotoxins.

Click [HERE](#) for ordering information and pictures of BASIK brand headspace vials.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024