



Are the BASIK brand headspace vials free of endotoxins – FAQ

These vials are not tested for the presence of endotoxins but since they are manufactured at 600°C, bacteria or organic matter will not likely exist after manufacturing.

While the vials are packed directly after cooling, they are *not* certified to be “untouched by human hands” and may or may not contain endotoxins.

Click [HERE](#) for ordering information and pictures of BASIK brand headspace vials.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com