



## Can the stainless steel tubing cutter wheel type accommodate different size tubing – FAQ

Yes, the stainless-steel tubing cutter (cat# 46001-01) has a screw which is tightened to secure the tubing and can be adjusted for a range of different tubing diameters. It can be used for any tubing diameter in the range of 1/32 inch (0.79 mm) to 1/4 inch (6.4 mm).



Click [HERE](#) for ordering information for the tubing cutter and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 03-05-2024