



How can I use the MicroSolv Window Maker to burn a larger window – FAQ

To begin, the window sizes are approximate and we cannot put a specification on exact widths. Since the user must hold the capillary in place and if it is not held steady, the window length may be longer.

Because of the approximate nature of the window length, which is generally sufficient, for users that want to make a larger **CE** window for which no size replacement module is available, it is advisable to use two or more portions of a smaller size to cover a large enough length of capillary.

Click [HERE](#) for Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com