



## What makes glass borosilicate glass – Tech Information

The glass material must contain at least 5% boron oxide to be termed “borosilicate glass”. Other types of glass with less than 5% is often called flint glass or sodium glass.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-03-2024