



Can the Window Maker heating modules be used with either US or European voltage units – FAQ

The Window Maker units are currently produced with 110/120V voltage for use in US & Canada (cat# [07200-S](#)) and the replacement modules (cat# [07301-S](#), [07201-S](#), [07204-S](#), [07207-S](#)) differ only in size of the burn window and can be used with new and older Window Makers even if you have a discontinued unit with different voltages.

Click [HERE](#) for MicroSolv Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com