

## Using an Aqueous Normal Phase ANP method will sodium ions interfere with succinate when analyzing sodium succinate - FAQ

***No they will not.***

Sodium ions, Na<sup>+</sup> will normally elute at or near the void volume in an Aqueous Normal Phase ANP method when analyzing polar compounds. Therefore, sodium should not interfere with quantitation of the succinate peak or any other peak for that matter.



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024