
White Paper: Effect of Diluent on Analyte Adsorption to RSA Glass Autosampler Vials.

It is well known that Silanol Groups present on the surface of Type I Glass Autosampler Vials contribute to quantitative loss in the case of Basic analytes. This is believed to occur primarily by electrostatic interaction and possibly by hydrogen bonding to some degree.

In order to investigate the nature of this interaction in more detail, we conducted a series of Experiments using different sample diluents.

Download the attached PDF to read full article using the link below.



[RSA Glass Vials Product Page](#)

Attachments

White Paper RSA Effect of Diluent on Adsorption Jan 14.pdf 2.3 Mb [Download File](#)