

## What is the chemical structure of the ligand in the Cogent DIOL HPLC stationary phase – FAQ

**The Cogent Diol™ stationary phase structure has a single point direct silicon to carbon attachment to the silica as shown below.**



**Note** this bonded phase is attached to the silica hydride surface and offers very unique retention mechanisms for both RP and ANP HPLC.



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-08-2024