

For short term reuse, store Simplus™ capillaries at room temperature. Many users simply tape them to the wall near the instrument.

For long term storage and reuse, fill with water and store at room temperature. Before use, store at room temperature. Never refrigerate or heat these before use or during storage.

Click [HERE](#) for capillaries ordering information

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)