

## Which HPLC column ID dimensions are the Hichrom brand guard column holders cat no 82000-GH compatible with - FAQ

The Hichrom brand guard column holder with the catalog number [82000-GH](#) can be used with any guard column in the range of 1.0mm – 4.6mm internal diameter (ID). We do not suggest using these guard columns with analytical columns outside of this size range.

In addition to the holder and guard column, you will need a coupler which is used to connect the guard column holder to the analytical column. [Click HERE](#) for ordering information and a picture of the union / column coupler.

For smaller ID columns we recommend using a Cogent Column Filter instead such as 49650-005. [Click HERE](#) for ordering information for 2um Cogent Column Filters.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)