



What are the dimensions of 4.6mm ID 5um frits for columns made with PEEK – FAQ

The Frit itself has an outer diameter of 0.188" and a thickness of 0.062" and is in a PEEK ring with 0.25" O.D. and 0.181" ID. The Frit is press-fit into the ring without glues or binders.

Click [HERE](#) for PEEK column hardware ordering page.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 10-22-2024