



## What are the dimensions of 4.6mm ID 5um frits for columns made with PEEK – FAQ

The Frit itself has an outer diameter of 0.188" and a thickness of 0.062" and is in a PEEK ring with 0.25" O.D. and 0.181" ID. The Frit is press-fit into the ring without glues or binders.

Click [HERE](#) for PEEK column hardware ordering page.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-05-2024