



Sample loop information - FAQ

Sample Loops are available in a variety of internal volumes.

Some of these will have several turns of the loop but others will not have a loop at all. In the latter case, the reason is because the tubing length is too small to make a loop. Whether or not there is a loop will not affect the functionality of the product.



Click [HERE](#) for sample loops ordering information and pictures.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com