

While BASIK™ brand screw caps are not designed for ion chromatography, they may work for some applications. But I would not use them sensitive methods.

We do suggest it might be worthwhile for maximum performance in IC to use our ultra-pure polyethylene screw caps, see image below.



Note that these caps, are a single use cap and will not reseal after injection.

Click [HERE](#) for more information about IC compatible caps.

Click [HERE](#) for PP vials that are IC compatible.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com