



## Do you have a tear off crimp cap that is also magnetic – FAQ

**No, we do not offer a magnetic tear off cap at this time because it requires aluminum (a soft metal) to allow easy tearing of the cap and aluminum is not magnetic.** The magnetic material used in crimp caps is not soft enough to tear.

Click [HERE](#) for crimp cap ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024