



## What is the maximum pressure tolerance for 0.020 ID stainless steel tubing - FAQ

**The pressure tolerance for this tubing is 15,800 psi.**

Please refer to our [maximum pressure tolerance chart](#) for full listings of our tubing products. Click [HERE](#) for high pressure tubing ordering information.

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)