



What is the maximum pressure tolerance for 0.020 ID stainless steel tubing – FAQ

The pressure tolerance for this tubing is 15,800 psi.

Please refer to our [maximum pressure tolerance chart](#) for full listings of our tubing products. Click [HERE](#) for high pressure tubing ordering information.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 10-11-2024