



## Which guard column do I use with the Partisil SAX Columns – FAQ

**We suggest you use catalog # 66900-AGC (Anion Guard Cartridge (5 pack)), for use with Partisil™ HPLC columns.**

You will also need a holder which is our catalog # 66900-WCS ,WCS guard cartridge holder, for use with Partisil HPLC columns.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-14-2024