



Which guard column do I use with the Partisil SAX Columns - FAQ

We suggest you use catalog # 66900-AGC (Anion Guard Cartridge (5 pack)), for use with Partisil™ HPLC columns.

You will also need a holder which is our catalog # 66900-WCS ,WCS guard cartridge holder, for use with Partisil HPLC columns.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024